



## Instructions

- 1 Print or trace mask onto Tyvek™ and cut center hole
- 2 Pre-fold mask on gray lines where indicated
- 3 Cut out shapes on solid black lines
- 4 Trace smaller pieces onto materials and cut
- 5 Peel & stick adhesive and nosepiece onto mask
- 6 Peel off paper backing and fold mask
- 7 Peel and stick seal along inside edge
- 8 Try on mask and mold nosipiece
- 9 Add strips of extra neoprene where there are gaps
- 10 Pop in 3D printed filter enclosure and add filter

